

# Complementa tu laboratorio con nuestro **ANALIZADOR VETERINARIO DE QUÍMICA SANGUÍNEA**



- Proceso de prueba sencillo.
- Resultados en 12 minutos.
- Control de calidad interno.
- Costo 0 de mantenimiento.



**DESEGO Panamá**

T: +5072898480/8481 C: +50763943784 [pty@desego.com](mailto:pty@desego.com)  
Avenida Balboa edicio H2O local D, Panamá.

**DESEGO**

# Veterinary Rotor Profile & Order Code

Conventional Test  
AW00245

Health Check  
AW00085

Liver Function  
AW00087

Kidney Function  
AW00198

Electrolyte  
AW00234

Primary Diagnostic  
AW00468

Pre-operation  
AW00235

Analyte	Full Name	Conventional Test AW00245	Health Check AW00085	Liver Function AW00087	Kidney Function AW00198	Electrolyte AW00234	Primary Diagnostic AW00468	Pre-operation AW00235
ALB	Albumin	●	●	●	●			
ALP	Alkaline phosphatase	●		●			●	●
ALT	Alanine aminotransferase	●	●	●			●	●
AMY	Amylase	●	●					
AST	Aspartate aminotransferase		●	●				●
Ca	Calcium		●		●	●		
CHE	Cholinesterase	●						
Cl	Chlorine					●		
CK	Creatine Kinase		●					●
CREA	Creatinine	●	●		●		●	●
GLU	Glucose	●	●		●		●	●
GGT	Gamma-glutamyltransferase			●				
GLOB*		●	●					
K	Potassium	●				●		
LDH	Lactate dehydrogenase							●
LPS	Lipase							
Mg	Magnesium					●		
Na	Sodium	●				●		
PHOS	Phosphorus		●		●	●		
TBA	Bile Acids			●				
TB	Total Bilirubin	●	●	●				
TC	Total Cholesterol			●				
Co2	Bicarbonate				●	●		
TG	Triglycerides		●					
TP	Total protein	●	●	●			●	●
UA	Uric Acid	●			●			
BUN	Urea nitrogen	●	●		●		●	●
CRP	C-reactive protein							

## The Special Biochemical Reagent Disc

Canine Inflammation  
AW00470

Coagulation  
AW00650

Analyte	Full Name	Canine Inflammation AW00470	Coagulation AW00650
LPS	lipase	●	
CRP	C-Reactive Protein	●	
AMY	Amylase	●	
UREA	Urea nitrogen	●	
CREA	Creatinine	●	
PT	Prothrombin		●
APTT	Activated partial thromboplastin		●
TT	Thrombin		●
FIB	Fibrinogen		●